

# THE ROCKER BUILT TO LAST

## THE INFINITY ROCKER™

THE LAST ROCKER YOU WILL EVER NEED TO BUY



~~1 Year  
Warranty~~

**4 YEAR  
WARRANTY**

This quiet, rugged laboratory rocker is ideal for all mixing applications. Have a shared rocker that's on almost all the time? The **INFINITY ROCKER™** is ideal for 24/7 operation as well as for occasional use.

### ■ 4 year complete warranty on motor and electronics

We've designed the **INFINITY ROCKER™** with constant, heavy-duty use in mind, and designed it to last - unlike all those rockers with one-year warranties.

### ■ Digital display of rocking speed

The **INFINITY ROCKER™** has a bright digital display that shows the rocking speed. Rocking range spans from 1 to 80 cycles per minute.

### ■ Available in multiple platform sizes and single, double, triple, or quadruple platform configurations

**INFINITY ROCKERS™** come with 11 x 14, 8 x 9, or 20 x 16 inch platforms. Stacking trays, Stack 3 and Stack 4 models provide extra capacity while keeping the footprint small.



Visit [www.nextadvance.com](http://www.nextadvance.com) or contact [info@nextadvance.com](mailto:info@nextadvance.com)

## NEXT >>> ADVANCE

INNOVATIVE PRODUCTS • BETTER WORKFLOW  
ACCURATE RESULTS • [www.nextadvance.com](http://www.nextadvance.com)

24 Prospect Avenue  
Averill Park, NY 12018, USA  
Telephone: 800-738-1681  
Fax: 518-674-0189  
[info@nextadvance.com](mailto:info@nextadvance.com)

# THE ROCKER BUILT TO LAST

## THE INFINITY ROCKER™

THE LAST ROCKER YOU WILL EVER NEED TO BUY

Available in five different models

### Infinity Rocker™ Pro

The basic rocker - rugged, reliable and quiet. 11 X 14 inch platform, single stack.  
Pro Stacking Tray available for double capacity.

### Infinity Rocker™ XL

Single stack rocker with a larger platform- 16 X 12 inch. XL Stacking Tray available for double capacity.

### Infinity Rocker™ Mini

Save space with the smallest model- 8 X 9 inch single platform. Mini Stacking Tray available for double capacity.

### Infinity Rocker™ Stack-3

Three times the capacity of the basic model. 11 X 14 inch platforms.  
Over 6.5 inches clearance between levels.

### Infinity Rocker™ Stack-4

Four times the capacity of the basic model. 11 X 14 inch platforms.  
Over 4 inches clearance between levels.

## Additional Features

- Adjustable tilt angle from 1° to 18°
- Quiet belt drive operation
- Improved rocking range spanning 1- 80 cycles per minute
- Safe for cold room or incubator use (4°C - 40°C)
- Microprocessor controlled motor speed
- Made in the USA

Visit [www.nextadvance.com](http://www.nextadvance.com) or contact [info@nextadvance.com](mailto:info@nextadvance.com)



**Pro Stacking Tray for Infinity Rocker™ Pro** - Snaps into place to double the capacity of your rocker

## NEXT >>> ADVANCE

INNOVATIVE PRODUCTS • BETTER WORKFLOW  
ACCURATE RESULTS • [www.nextadvance.com](http://www.nextadvance.com)

24 Prospect Avenue  
Averill Park, NY 12018, USA  
Telephone: 800-738-1681  
Fax: 518-674-0189  
[info@nextadvance.com](mailto:info@nextadvance.com)