

# YOUR NEXT ROCKER

## THE NEXT ROCKER



- > QUIET
- > RUGGED
- > RISK FREE

### > Adjustable Tilt

The tilt angle adjustment screw in the back of each rocker easily lets you change the maximum tilt angle from 1° to ± 18°

### > Stacking Trays

Double the capacity with TIER 2 stacking trays. They simply snap onto the platform; their spring loaded legs hold them securely in place.

### > Choose your platform size

Pro: 11 x 14 inch (27 x 35 cm)  
Mini: 8 x 9 inch (20 x 23 cm)

### > High Value

The simple reliable designs and Next Advance's efficiency enable these rockers to sell for much less than comparable models from other companies.

## NEXT >>> ADVANCE

INNOVATIVE PRODUCTS • BETTER WORKFLOW  
ACCURATE RESULTS • [www.nextadvance.com](http://www.nextadvance.com)

24 Prospect Avenue  
Averill Park, NY 12018, USA  
Telephone: 800-738-1681  
Fax: 518-674-0189  
[info@nextadvance.com](mailto:info@nextadvance.com)

# YOUR NEXT ROCKER

The **NEXT ROCKER PRO** (right) is the standard NEXT ROCKER model with a 11 x 14 inch platform. The NEXT ROCKER MINI (8 x 10 inches) is small enough to fit on a sheet of paper.



You can place a stacking tray on a platform (below) or remove it when desired without any tools. The spring-loaded legs hold the stacking tray in place.



**TIER 2 Stacking Trays** (left) secure onto the platforms of the NEXT ROCKERS. This gives you the ability to expand your rocker 4.5 inches above the main platform. Easy to clean and remove!



## > RUGGED CONSTRUCTION

These sturdily built rockers sport a stainless steel front panel and platform and a heavy gauge enclosure. A top grade Pittman instrument motor drives the rocking and a pulse width modulation circuit controls the rocking speed.

## > QUIET OPERATION

NEXT ROCKERS are very quiet due to the high quality motor and vibration isolation mounts.

## > RISK FREE PURCHASE

Every NEXT ROCKER comes with a 30 day money back guarantee.

## > FEATURES

- Speed ranges up to 80 cycles/minute (160 rock actions per minute)
- Tilt angle adjustable from 1° to ± 18°
- Ideal for all mixing applications
- 11 x 14 inch (28 x 36 cm) stainless steel platform with a rubber non-slip pad
- Smaller (8 x 10 inch) rocker platform available
- Operational in constant 4° to 60°C environment
- CE/UL certified 24 volt power supply with plug configuration for wall outlets worldwide
- 1 year warranty
- 30 day money back guarantee

# NEXT >>> ADVANCE

INNOVATIVE PRODUCTS • BETTER WORKFLOW  
ACCURATE RESULTS • [www.nextadvance.com](http://www.nextadvance.com)

24 Prospect Avenue  
Averill Park, NY 12018, USA  
Telephone: 800-738-1681  
Fax: 518-674-0189  
[info@nextadvance.com](mailto:info@nextadvance.com)