



# 10cm<sup>2</sup> Tiny TISSUE CULTURE FLASKS



- Total Access to Cells
- Ideal for Primary Cell Lines, Organ Harvesting, Whole Tissue Growth, Colony Picking, etc!
- Limit Transfer Steps

## Key Points

Outstanding for Small or Primary Cultures



Ideal for Cell Scraping, Colony Picking, or Cell Lifting



Stand for Stability on Shakers



The "Tiny" Flask Viewing Area



The "Tiny" Flask eliminates the Transfer Step; It Fits Most 50ml Rotors



| Product # | Version                         | Growth surface     | Lid opening | Volume     | Dimensions  | Pieces / Pack | Pieces / Case | Your Price |
|-----------|---------------------------------|--------------------|-------------|------------|-------------|---------------|---------------|------------|
| 91253     | 10 cm <sup>2</sup> "Tiny" Flask | 10 cm <sup>2</sup> | 2.6 cm      | Max of 6ml | 30 x 110 mm | 5             | 350           | \$229.00   |

Available from:

**NEXT » ADVANCE**

Innovative Lab Products for the Life Sciences

24 Prospect Ave  
Averill Park, NY 12018

(800) 738-1681  
[www.nextadvance.com](http://www.nextadvance.com)