

# MAKE CONNECTIONS!

## ONE>CLICK™ BOTTLE CAPS



The easiest fluidic connections to laboratory bottles and carboys.

- > QUICK CONNECT & DISCONNECT
- > FREE TO SWIVEL, AVOIDS KINKS
- > LEAK-FREE

### > QUICK CONNECT & DISCONNECT


Just click the tubing onto the bottle cap to connect. To disconnect, just click the tab on the fitting and slip the tubing off the cap. No more tugging.

### > FREE TO SWIVEL, AVOIDS KINKS

The tubing is free to pivot so the tubing will not twist or coil. Bottle caps can be screwed on or off without lifting the bottle or twisting the tubing.

### > LEAK-FREE

An O-ring seal prevents leakage at the connection. In addition, a valve inside the cap prevents fluids from escaping when tubing is not connected.



ONE>CLICK™ BOTTLE CAPS ARE AVAILABLE FOR MANY LABORATORY BOTTLES AND CARBOYS.  
CAPS ARE SUPPLIED WITH BARBED FITTINGS THAT SLIP INTO THE END OF THE TUBING.

## NEXT>>>ADVANCE

24 Prospect Avenue, Averill Park, NY 12018, USA | Telephone: 1.800.738.1681 | Fax: 1.610.862.1718  
email: [info@nextadvance.com](mailto:info@nextadvance.com) | [www.nextadvance.com](http://www.nextadvance.com)