


*Disrupt, homogenize,  
dissociate or lyse cells,  
tissues and plant matter*



## Stomp your samples flat In seconds. Right in the tube.

### FEATURES

- ▶ Can accommodate 5 mL and microcentrifuge tubes.
- ▶ Powerful lever action makes even the toughest samples easy to crush.
- ▶ Quick and easy. Stomping only takes a few seconds.
- ▶ Disposable pestle shield prevents cross-contamination.



### Your problem samples won't know what hit them

With the Next Advance Stomper™, you can crush tough material in seconds, with no cross-contamination.

Use it with the Bullet Blender® homogenizer to process samples like soybeans, shells, and corn!

#### Collect



#### Stomp



#### Homogenize



**Step 1:** Place your sample in the tube.

**Step 2:** Stomp it.

**Step 3:** Add grinding beads and homogenize in the Bullet Blender®. Your sample will be a fine powder in as *little as 5 minutes*.

### For More Information, Contact:

**Next Advance, Inc.**

Tel. 1.800.738.1681

Fax. 1.518.674.0189

info@nextadvance.com

www.nextadvance.com