

Freedom Rocker™ BlotBot® 240 Specifications

Capacity:

- Up to 2 Midi Blots or 4 Mini Blots.
- Up to 2 reagents per blot, up to 50 mL per reagent per blot, not including buffer
- Trays hold up to 75 mL of buffer per miniblot, 150 mL per midiblot.
- Includes 2L bottle for buffer, option for up to 4L

Features:

- Minimum of 3.5 mL of antibody solution per miniblot
- All reagents can be recovered.
- Stores up to 8 protocols
- Protocols can include up to 9 steps, with up to 10 washes per step
- Each step can include different volumes and rocking speeds
- Rocking speeds adjustable from 5 to 40 cycles per minute.
- Automatic restart after power outage
- Time per incubation or wash adjustable from 1 minute to 96 hours

General:

- Dimensions: 60 cm wide by 41 cm deep by 45 cm high.
- Relative Humidity: 5- 90% non-condensing.
- Power Requirement: 12 VDC, 5 Amp.
- Can be used in temperatures ranging from 4°C to 40°C.
- Altitude: <2000m..
- Storage Temperature: -40°C to 50°C
- 2 year warranty.
- Made in the USA.