

## ZROB015

Zirconium oxide grinding media provided by Next Advance, Inc. is ideal for the homogenization of biological samples. Due to the lower density of the material, it is suitable for the homogenization of softer samples, such as brain or kidney tissue. Available as RNase-free or non-sterile.

RNase-free zirconium oxide beads (ZROB015-RNA) undergo a special cleaning process. They are guaranteed and certified to be RNase-free.



### Uses

- Homogenize tissue to release protein, DNA, RNA or other molecules of interest.
- Release viable bacteria from tissue.
- Release viable cells from tissue to generate suspensions.

### Properties

Description	Property
Shape and Appearance	Spherical and white in color
Composition	95% Zirconium Oxide & 5% Yttrium Oxide
Specific Density	Greater than 5 g/cm <sup>3</sup>
Diameter	0.15 mm ± 0.02 mm