

Overview

Zirconium oxide grinding media provided by Next Advance is ideal for the homogenization of biological samples. Due to the lower density of the material, it is suitable for the homogenization of softer samples, such as brain or kidney tissue.

Uses

- Homogenize tissue to release protein, DNA, RNA or other molecules of interest.
- Release viable cells from tissue to generate suspensions.

Properties



Description	Property
Shape	Spherical
Composition	95% Zirconium Oxide 5% Yttrium Oxide
Specific Density	~ 6.0 g/cm ³
Bulk Density	~ 3.75 g/cm ³
Diameter	6.0 mm ± 0.10 mm

MK19-0261A