

Glass Beads, 0.1 mm

Overview

Glass grinding media provided by Next Advance is ideal for the homogenization of biological samples. The low density of the material makes it suitable for the homogenization of soft tissues and easy-to-lyse cell culture.

Uses

- General-purpose bead beating in bead mill homogenizers for small soft biological samples.
- Homogenize samples to release DNA or other molecules of interest.
- Compatible with all Bullet Blender homogenization models.

Example Sample Types

- Bacteria and small cell culture.
- Cells with thin or no cell walls.
- Including other samples with similar characteristics.



Specifications

Product Amount	1 lb.
Shape	Spherical
Material	Glass
Composition	Type S lead free glass
Specific Density	~ 2.5 g/cm ³
Bulk Density	~ 1.5 g/cm ³
Average Size	0.1 mm ± 0.02 mm

Usage Statement

This product is intended solely for research and laboratory use. It is not intended for use in humans or animals, nor for any clinical, diagnostic, or therapeutic applications. Users are responsible for ensuring compliance with all applicable standards and regulatory requirements relevant to their specific workflows.

MK19-0242A