

Stainless Steel Beads, 5.6 mm

Overview

Stainless steel grinding media provided by Next Advance is ideal for the homogenization of biological samples. The unique UFO shape and higher density of the material make it suitable for the homogenization of tough samples.

Uses

- General-purpose bead beating in bead mill homogenizers for tough samples.
- Powderizing dry samples.
- Homogenizing samples to release bacteria or organelles from cells.
- Compatible with 5E or 50DX Bullet Blender homogenization models.



Example Sample Types

- Tough tissue like heart or skin samples.
- Fibrous plant matter like leaf or root.
- Small insects like tick or mosquito.
- Soil or fecal matter.
- Including other samples with similar characteristics.

Specifications

Product Amount	50 mL
Shape	UFO shaped
Material	Stainless Steel
Composition	300 series stainless steel
Specific Density	~ 7.8 g/cm ³
Bulk Density	~ 5.5 g/cm ³
Average Size	5.6 mm ± 0.08 mm

Usage Statement

This product is intended solely for research and laboratory use. It is not intended for use in humans or animals, nor for any clinical, diagnostic, or therapeutic applications. Users are responsible for ensuring compliance with all applicable standards and regulatory requirements relevant to their specific workflows.