

Zirconium Oxide Beads, 2.5 mm

Overview

Zirconium oxide grinding media provided by Next Advance is ideal for the homogenization of biological samples. The size and density of the material make it suitable for the homogenization of tissue with soft-medium toughness and for releasing bacteria from tissue.

Uses

- General-purpose bead beating in bead mill homogenizers for tissues with soft-medium toughness.
- Homogenize samples to release bacteria or organelles from cells.
- Homogenize samples to generate cell suspensions.
- Compatible with all Bullet Blender homogenization models.



Example Sample Types

- Tissues with soft-medium toughness like liver, brain, or kidney.
- Including other samples with similar characteristics.

Specifications

| | |
|-------------------------|--------------------------------------|
| Product Amount | 1 lb. |
| Shape | Spherical |
| Material | Zirconium Oxide |
| Composition | 95% Zirconium Oxide 5% Yttrium Oxide |
| Specific Density | ~ 6.0 g/cm ³ |
| Bulk Density | ~ 3.75 g/cm ³ |
| Average Size | 2.5 mm ± 0.08 mm |

Usage Statement

This product is intended solely for research and laboratory use. It is not intended for use in humans or animals, nor for any clinical, diagnostic, or therapeutic applications. Users are responsible for ensuring compliance with all applicable standards and regulatory requirements relevant to their specific workflows.